

# Quality Scores for Next-Generation Sequencing

Assessing sequencing accuracy using Phred quality scoring.

## Introduction

A next-generation sequencing experiment consists of a series of discrete steps that uniquely contribute to the overall quality of a data set. Sequencing quality metrics can provide important information about the accuracy of each step in this process, including library preparation, base calling, read alignment, and variant calling. Base calling accuracy, measured by the Phred quality score (Q score), is the most common metric used to assess the accuracy of a sequencing platform. It indicates the probability that a given base is called incorrectly by the sequencer.

Historically used to determine Sanger sequencing accuracy, Phred originated as an algorithmic approach that considered Sanger sequencing metrics, such as peak resolution and shape, and linked them to known sequence accuracy through large multivariate lookup tables. This method proved to be highly accurate<sup>1</sup> across a range of sequencing chemistries and instruments, making it the quality scoring standard for commercial sequencing technologies.

While next-generation sequencing metrics vary from those of Sanger sequencing (e.g., no electropherogram peak heights), the process of generating a Phred quality scoring scheme is largely the same. Parameters relevant to a particular sequencing chemistry are analyzed for a large empirical data set of known accuracy. The resulting quality score lookup tables are used to calculate a quality score for *de novo* next-generation sequencing data (in real time on Illumina platforms), possessing an equivalent meaning to the historical metrics familiar to most Sanger sequencing users.

## Calculating Phred Quality Scores

Q scores are defined as a property that is logarithmically related to the base calling error probabilities ( $P$ )<sup>2</sup>.


$$Q = - 10 \log_{10} P$$

For example, if Phred assigns a Q score of 30 (Q30) to a base, this is equivalent to the probability of an incorrect base call 1 in 1000 times (Table 1). This means that the base call accuracy (i.e., the probability of a correct base call) is 99.9%. A lower base call accuracy of 99% (Q20) will have an incorrect base call probability of 1 in 100, meaning that every 100 bp sequencing read will likely contain an error. When sequencing quality reaches Q30, virtually all of the reads will be perfect, having zero errors and ambiguities. This is why Q30 is considered a benchmark for quality in next-generation sequencing. By comparison, Sanger sequencing systems generally produce base call accuracy of ~99.4%, or ~Q20<sup>3</sup>. Low Q scores can increase false-positive variant calls, which can result in inaccurate conclusions and higher costs for validation experiments.

## Illumina Data Quality

Illumina Q score calculations have been shown to be very similar to the actual data quality observed in human genome sequencing<sup>4</sup>. Figure 1 shows that predicted and empirical quality scores from a HiSeq 2000

Figure 1: High Correlation of Empirical and Predicted Q Scores



Illumina sequencing Q scores are highly accurate. This example shows that predicted Q scores for a HiSeq 2000 run correlate well to empirically derived Q scores.

run are well correlated. Q scores can reveal how much of the data from a given run is usable in a resequencing or assembly experiment. Sequencing data with lower quality scores can result in a significant portion of the reads being unusable, resulting in wasted time and expense. PhiX quality scores for the MiSeq<sup>®</sup> and HiSeq<sup>®</sup> systems show that nearly all bases have scores > Q30 for single and paired-end reads (Figure 2). Comparison of *E. coli* whole-genome sequencing data shows that this high data quality is consistent across both platforms (Table 2).

Table 1: Quality Scores and Base Calling Accuracy

Phred Quality Score	Probability of Incorrect Base Call	Base Call Accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1,000	99.9%
40	1 in 10,000	99.99%
50	1 in 100,000	99.999%

Table 2: MiSeq vs HiSeq 2000 *E.coli* K12 MG1655 Data Comparison

Metric	MiSeq System		HiSeq System	
	Read 1	Read 2	Read 1	Read 2
% Bases Q ≥ 30	91.9	87.5	89.3	86.1
% Total Bases Q ≥ 30	89.7		87.7	

A whole-genome sequencing run (2 × 150 bp) of *E. coli* K12 MG1655 performed on the MiSeq system yielded 1.7 Gb of high-quality data. MiSeq data were trimmed to 2 × 100 bp to allow for a direct comparison with 2 × 100 bp reads from the HiSeq 2000 platform.

